



wonderbrush

multi-purpose soldering brush!

- * For use after laser, fazer, and torch soldering to wipe off the oxides or the flux
- * For use on all metals, including precious, non- precious, and even titanium
- * For use on partials
- * For use on chrome cobalt



dentalcreations ltd.

www.wonderfill.com | +1.254.772.4661

DENTAL CREATIONS, LTD., P.O. BOX 21325, WACO, TEXAS 76702 U.S.A.

DIRECTIONS FOR USE:

After soldering, to wipe off the oxides rub the area with the soldering brush. Before applying to porcelain, it is recommended to lightly blast with aluminum oxide.