

Products

WonderFill

- Introduction
- The Story
- How it Works
- Compare
- Free Sample

WonderForm

- Product Details
- How it Works

WonderOrtho

- [Product Details](#)

WonderAdmix

- Product Details
- Free Sample

WonderPeg

- Product Details
- How it Works
- Use It For . . .
- Free Sample

wonderortho

Perfecting the process of making study models and it's EZ.
Produce better study models in less time.

The box and pour process has never been so easy.

The WonderOrtho simplified application eliminates all mixing of materials, eliminates all wax products, and dramatically reduces the time to prepare cast models. Research shows that using the innovative WonderOrtho process saves an average of 25 minutes per model.



BENEFITS INCLUDE:

**No Mixing, No Mess,
No Wax Products**

Works with ALL Laboratory, Orthodontic, and General Dentistry Applications

Works with ANY Casting Material, Stone, Silicone, and Rubber-Based products

Results in a Superior Quality Mold while Taking Less than Five Minutes to Complete

254.772.4661